


Product Datasheet

SA2 antibody (orb704237)

Description	SA2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human SA2 (1021-1120/1231 aa)
Target	SA2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	141 kDa
Uniprot ID	Q8N3U4
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: NIH/3T3(Mouse)
Cell;...