


Product Datasheet

CD34 antibody (orb704258)

Description	CD34 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human CD34
Target	CD34
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	39 kDa
Uniprot ID	P28906
Dilution Range	IHC-P=1:50-200, IHC-F=1:50-200, ICC=1:50-200, IF=1:50-200, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse testis
using CD3...