



Product Datasheet

CYP26A1 antibody (orb704276)

Description	CYP26A1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human CYP26A1 (200-400/497 aa)
Target	CYP26A1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	48 kDa
Uniprot ID	O43174
Dilution Range	WB=1:500-1000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human liver
carcinoma ...



Western
blot: Lane
1: HepG2
(Human)
Cell; L...