


Product Datasheet

C Met antibody (orb704283)

Description	C Met antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human MET (620-720/1390 aa)
Target	MET
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33/123/153 kDa
Uniprot ID	P08581
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Mouse Kidney; Lane 2: ...