



Product Datasheet

BST2 antibody (orb704295)

Description	BST2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human BST2
Target	BST2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	18 kDa
Uniprot ID	Q10589
Dilution Range	WB=1:300-500, IHC-P=1:50-200, Flow-Cyt=1:50, IHC-F=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human colon
cancer usi...



Western
blot:Lane
1:Hela
Cell;Lane
2:Sih...