

Product Datasheet

ARP3 antibody (orb704316)

Description	ARP3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant protein of C-terminal human Arp3.
Target	ARP3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	47 kDa
Uniprot ID	P61158
Dilution Range	WB=1:500-1000, IHC-P=1:50-200, ICC=1:50-100, Flow-Cyt=1:50-100, IF=1:50-100
Expiration Date	12 months from date of receipt.