



Product Datasheet

MKLP1 antibody (orb704336)


Description	MKLP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human MKLP1
Target	MKLP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	110 kDa
Uniprot ID	Q02241
Dilution Range	WB=1:500-1000, ICC=1:50, Flow-Cyt=2ug/Test, IF=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometry of A549. Primary Antibody ...



Western blot: Lane 1:293 Cell; Lane 2:A549...



Western blot: Lane 1:Rat Lymph node; Lane ...