


Product Datasheet

Collagen I antibody (orb704339)

Description	Collagen I antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Collagen I
Target	Collagen I
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	130 kDa
Uniprot ID	P02452
Dilution Range	IHC-P=1:50-200, IF=1:100-500, WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
kidney tissue
us...