


Product Datasheet

Cleaved PARP antibody (orb704341)

Description	Cleaved PARP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant human PARP protein
Target	Cleaved PARP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	89/116 kDa
Uniprot ID	P09874
Dilution Range	WB=1:500-2000, Flow-Cyt=1:50, ICC=1:50-200
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Jurkat Cell; Lane 2: A...