



Product Datasheet

ACE2 antibody (orb704343)


Description	ACE2 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ACE2 (631-805/805 aa)
Target	ACE2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	87 kDa
Uniprot ID	Q9BYF1
Dilution Range	WB=1:500-1000, IHC-P=1:400-800, ICC=1:100
Expiration Date	12 months from date of receipt.



ICC analysis of SH-SY5Y cell using ACE2 ...



IHC-P analysis of Human kidney using ACE...



Western blot: Lane 1: Recombinant human AC...