

Product Datasheet

LYVE1 antibody (orb704345)

Description	LYVE1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	Recombinant human LYVE1 protein
Target	LYVE1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	32 kDa
Uniprot ID	Q9Y5Y7
Dilution Range	IHC-P=1:50-200, ICC=1:50, Flow-Cyt=1:50, IHC-F=1:50-200
Expiration Date	12 months from date of receipt.