

Product Datasheet

phospho C-Myc (Thr58 antibody (orb704368))

Description	phospho C-Myc (Thr58 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Myc around the phosphorylation site of Thr58 and Ser62 (LP(p-T)PPL(p-S)PS)
Target	phospho C-Myc (Thr58 / Ser62)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	P01106
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.