



Product Datasheet

SQSTM1 antibody (orb704394)

Description	SQSTM1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human SQSTM1/p62 (1-440 aa)
Target	SQSTM1/p62
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38/47 kDa
CAS Number	Sequestosome-1
Uniprot ID	Q13501
Dilution Range	WB=1:500-2000, IHC-P=1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-SQSTM1/...



Western blot: Lane 1:293T (Human) Cell;La...