



Product Datasheet

YY1 antibody (orb704397)


Description	YY1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant human YY1
Target	YY1(Nuclear Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	46 kDa
Uniprot ID	P25490
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, Flow-Cyt=1ug/Test, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometry of K562.Primary Antibody ...



IHC-P analysis of Mouse bladder using YY...



IHC-P analysis of Human smooth muscle us...