


Product Datasheet

Caspase-9 antibody (orb704406)

Description	Caspase-9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Caspase-9
Target	Caspase-9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	35/50 kDa
Uniprot ID	P55211
Dilution Range	WB=1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Hela (Human) Cell; La...