


Product Datasheet

Bcl-2 Beta antibody (orb704452)

Description	Bcl-2 Beta antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human apoptosis regulator Bcl-2 beta isoform (51-150/205 aa)
Target	Bcl-2 beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	P49950
Dilution Range	Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow
cytometry of
HL-
60.Primary
Antibody...