



Product Datasheet

SNAIL + SLUG (phospho-Ser246) antibody (orb704478)

Description	SNAIL + SLUG (phospho-Ser246) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human SNAIL around the phosphorylation site of Ser246 (TF(p-S)RM)
Target	phospho-SNAIL + SLUG (Ser246)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	O95863
Dilution Range	IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
colon
carcinoma
...



IHC-P
analysis of
rat heart
using
phosph...