



Product Datasheet

FoxP3 antibody (orb704485)


Description	FoxP3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human FoxP3 (351-431/431 aa)
Target	FOXP3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	Q9BZS1
Dilution Range	WB=1:500-2000, ICC=1:100, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of MCF7.Primary Antibody ...



ICC analysis of U2OS cell using FoxP3 po...



Western blot:Lane 1:Thymus (Mouse);Lane ...