



Product Datasheet

Caspase-1 P10 antibody (orb704487)

Description	Caspase-1 P10 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse Caspase-1 P10 (301-380/402 aa)
Target	Caspase-1 P10
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10/45 kDa
Uniprot ID	P29466
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Sample:
Lane 1:
Small
intestine
(Rat) Ly...



Blank
control:
Raw264.7.
Primary
Antibod...