

Product Datasheet

NET1 antibody (orb704493)

Description	NET1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human NET1 (501-596/596 aa)
Target	NET1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	Q7Z628
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:25, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.