



Product Datasheet

BANF1 antibody (orb704500)


Description	BANF1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human BANF1 (Full length)
Target	BANF1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	10 kDa
Uniprot ID	075531
Dilution Range	IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
esophageal
cance...



IHC-P
analysis of
rat testis
using
BANF1...



IHC-P
analysis of
mouse testis
using BAN...