



Product Datasheet

CD38 antibody (orb704512)

Description	CD38 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human CD38
Target	CD38
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	34 kDa
Uniprot ID	P28907
Dilution Range	WB=1:500-1000, IHC-P=1:100-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
gastric
carcinom...



Western
blot:Primary
antibody at
1:2000 ...