



Product Datasheet

SLC15A1 antibody (orb704521)


Description	SLC15A1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SLC15A1
Target	SLC15A1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	79 kDa
Uniprot ID	P46059
Dilution Range	WB=1:500-1000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Liver (Mouse); Lane 2...



IHC-P analysis of rat kidney using SLC15...



IHC-P analysis of human colon carcinoma ...