



Product Datasheet

Integrin Alpha 2 antibody (orb704523)


Description	Integrin Alpha 2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Integrin alpha 2
Target	Integrin alpha 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	127 kDa
Uniprot ID	P17301
Dilution Range	WB=1:500-2000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of human gastric carcinom...



IHC-P analysis of rat colon using Integr...



Western blot analysis using Anti-Integri...