



Product Datasheet

ERK1 antibody (orb704524)


Description	ERK1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ERK1/2
Target	ERK1/2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	42/44 kDa
Uniprot ID	Q63844
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:100, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of Hela.Primary Antibody ...



ICC analysis of A549 cell using ERK1/2 m...



IHC-P analysis of rat colon using ERK1 2...