



Product Datasheet

RNF43 antibody (orb704525)


Description	RNF43 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human RNF43 (401-500/783 aa)
Target	RNF43
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	83 kDa
Uniprot ID	Q68DV7
Dilution Range	WB=1:500-2000, IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat brain
using RNF43
...



IHC-P
analysis of
human brain
using RNF4...



IHC-P
analysis of
mouse
pancreas
using R...