



Product Datasheet

HMGB1 antibody (orb704527)


Description	HMGB1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human HMGB1
Target	HMGB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	P09429
Dilution Range	WB=1:500-1000, IHC-P=1:200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat small
intestine us...



IHC-P
analysis of
human
gastric
carcinom...



IHC-P
analysis of
mouse lung
using
HMGB1...