



Product Datasheet

Alpha Smooth Muscle Actin antibody (orb704531)


Description	Alpha Smooth Muscle Actin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Actin alpha
Target	alpha smooth muscle Actin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	42 kDa
Uniprot ID	P62736
Dilution Range	WB=1:1000-5000, IHC-P=1:50-200, ICC=1:100
Expiration Date	12 months from date of receipt.



ICC analysis of HepG2 cell using alpha s...



IHC-P analysis of mouse uterus using alp...



IHC-P analysis of rat colon using alpha ...