



Product Datasheet

GAPDH (Loading Control) antibody (orb704534)


Description	GAPDH (Loading Control) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GAPDH (Loading Control)
Target	GAPDH (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	38 kDa
Uniprot ID	P04406
Dilution Range	WB=1:5000-50000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat testis
using
GAPDH...



IHC-P
analysis of
human liver
using
GAPD...



IHC-P
analysis of
mouse
testis using
GAP...