



Product Datasheet

PPM1D antibody (orb704535)


Description	PPM1D antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PPM1D (471-570/605 aa)
Target	PPM1D
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67 kDa
Uniprot ID	O15297
Dilution Range	IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse brain
using PPM1...



IHC-P
analysis of
rat kidney
using
PPM1D...



IHC-P
analysis of
human rectal
carcinoma...