

Product Datasheet

Mouse Gnai2 protein (orb704667)

Description

Recombinant Mouse Guanine nucleotide-binding protein G(i) subunit

Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Protein Sequence	GCTVSAEDKAAERSKMIDKNLREDGEKAAREVKLLLLGAGESGKSTIVKQMKII HEDGYSEEECRQYRAVVYSNTIQSIMAIVKAMGNLQIDFADPQRADDARQLFAL SCAAEEQGMLPEDLSGVIRRLWADHGVQACFGRSREYQLNDSAAYYLNDLERIA QSDYIPTQQDVLRTVKTGIVETHFTFKDLHFKMFDVGGQRSEKWKWICFEG VTAIIFCVALSAYDLVLAEDEEMNRMHESMKLFDSICNNKWFTDTSIILFLNKKDL FEEKITQSSLTICFPEYTGANKYDEAASYIQSKFEDLNKRKDTKEIYTHFTCATDTK NVQFVFDAVTDVVIKNNLKDCGLF
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Baculovirus
MW	44.2 kDa
Uniprot ID	P08752
Expiration Date	6 months from date of receipt.