

## Product Datasheet

### Human CCDC112 protein (orb704829)

## Description

Recombinant Human Coiled-coil domain-containing protein

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Endotoxins</b>        | Not test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Preservatives</b>     | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.                                                                                                                                                                                                                                                              |
| <b>Form/Appearance</b>   | Liquid or Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                 |
| <b>Tag</b>               | N-terminal 10xHis-tagged and C-terminal Myc-tagged                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Application notes</b> | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                     |
| <b>Protein Sequence</b>  | MEKDKHSHFYNQKSDFRIEHSMLEELLENKLIHSRTERAKIQQLAKIHNNVKKL<br>QHQLKDVKPTPDFVEKLREMMEEIENAINTFKEEQRLIYEELIKEEKTNNELSAIS<br>RKIDTWALGNSETEKAFRAISSKVPVDKVTSTLPEEVLD FEKFLQQTGGRQGA<br>WDDYDHQNFVKVRNKHKGKPTFMEEVLEHLPGKTQDEVQQHEKWYQKFLALE<br>ERKKESIQIWTKKQKQKREEIFKLKEKADNTPVLFHNKQEDNQKQKEEQRKKQKL<br>AVEAWKKQKSIEMSMKCASQLKEEEEKEKKHQKERQRQFKLKLLESYTQQKKE<br>QEEFLRLEKEIREKA EKA EKRNAADEISR FQERDLHKLELKILDRQAKED EKSQK<br>QRR LAKLKEKVENNVSRDPSRLYKPTKGWEERTKKIGPTGSGPLLHIPHRAIPTW<br>RQGIQRRV |
| <b>Purity</b>            | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>MW</b>                | 61.0 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Uniprot ID</b>        | <b>Q8NEF3</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |