



Product Datasheet

GAPDH antibody (orb705357)


Description	Mouse monoclonal antibody to GAPDH
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant protein corresponding to GAPDH.
Target	GAPDH
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 0.2% rAlbumin, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 0.2% rAlbumin, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P16858 , P04406 , P04797
Entrez	100042025 , 2597 , 24383
Dilution Range	WB: 1:3000-1:10000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of GAPDH expression...



Immunohistochemical analysis of GAPDH staining...



Immunofluorescent analysis of GAPDH staining...