



Product Datasheet

Topoisomerase 2 alpha antibody (orb705505)


Description	Mouse monoclonal antibody to Topoisomerase 2 alpha
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Topoisomerase 2 alpha. The exact sequence is proprietary.
Target	TOP2A
Preservatives	Mouse IgG2b. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Form/Appearance	Mouse IgG2b. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P11388
Entrez	7153
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:300, IF: 1:100-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of Topoisom...



Immunohistochemical analysis of Topoisom...



Immunohistochemical analysis of Topoisom...