



Product Datasheet

CD274 antibody (orb705620)


Description	Rabbit polyclonal antibody to CD274
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human CD274. The exact sequence is proprietary.
Target	CD274
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9NZQ7 , Q9EP73
Entrez	60533 , 29126
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of CD274 expression...



Immunohistochemical analysis of CD274 staining...



Immunofluorescent analysis of CD274 staining...