



Product Datasheet

SMAD5 antibody (orb705630)

Description	Rabbit polyclonal antibody to SMAD5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant fusion protein of human SMAD5. The exact sequence is proprietary.
Target	SMAD5
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q99717 , P97454 , Q9R1V3
Entrez	4090 , 59328 , 17129
Dilution Range	WB: 1:500-1:1000, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of SMAD5 expression...



Immunohistochemical analysis of SMAD5 staining...