

## Product Datasheet

### Rabbit Anti-Sheep IgG (H&L) antibody (Cy3) (orb705790)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal Secondary antibody against Rabbit Anti-Sheep IgG (H&L) conjugated to Cy3                           |
| <b>Species/Host</b>        | Rabbit                                                                                                               |
| <b>Reactivity</b>          | Sheep                                                                                                                |
| <b>Conjugation</b>         | Cy3                                                                                                                  |
| <b>Tested Applications</b> | ELISA, FC, IF, IH                                                                                                    |
| <b>Immunogen</b>           | Sheep IgG                                                                                                            |
| <b>Preservatives</b>       | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide          |
| <b>Form/Appearance</b>     | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                     |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                           |
| <b>Source</b>              | Rabbit                                                                                                               |
| <b>Dilution Range</b>      | IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                      |