

Product Datasheet

Rabbit Anti-Human IgG F(ab')₂ antibody (Cy5.5) (orb705889)

Description	Rabbit polyclonal Secondary antibody against Rabbit Anti-Human IgG F(ab') ₂ conjugated to Cy5.5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Cy5.5
Tested Applications	ELISA, FC, IF, IH
Immunogen	Human IgG F(ab') ₂
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
Expiration Date	12 months from date of receipt.