

Product Datasheet

Myc-tag antibody (FITC) (orb705936)

Description	Rabbit polyclonal antibody against Anti-Myc-tag conjugated to FITC
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of Myc-tag.
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	WB: 1:2000-1:5000, IF: 1:500-1:2000, FC: 1:1000-1:4000
Expiration Date	12 months from date of receipt.