

Product Datasheet

PKC beta (phospho-T641) Blocking Peptide (orb706051)

Description The peptide is used to block Anti-PKC beta (pT641) antibody <a href="https://www.biorbyt.com/prkcb-p..."

Reactivity Human, Mouse, Rat, Zebrafish

Tested Applications Blocking

Target PRKCB

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **P05771, P68404**

Entrez **18751, 5579**

Expiration Date 12 months from date of receipt.