

Product Datasheet

INT2 Blocking Peptide (orb706095)

Description The peptide is used to block Anti-INT2 antibody
<a href="https://www.biorbyt.com/ints2-antibody-orb5..."

Reactivity Human, Mouse, Rat

Tested Applications Blocking

Target INTS2

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **Q9H0H0, Q80UK8**

Entrez **57508, 70422**

Expiration Date 12 months from date of receipt.