

Product Datasheet

Serpin B12 Blocking Peptide (orb706143)

Description	The peptide is used to block Anti-Serpin B12 antibody <a href="https://www.biorbyt.com/serpinb12-ant..."
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	SERPINB12
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9D7P9, Q96P63
Entrez	89777, 71869
Expiration Date	12 months from date of receipt.