

Product Datasheet

PKC-pan Blocking Peptide (orb706237)

Description	The peptide is used to block Anti-PKC-pan antibody
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	PRKCA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q02156, P20444, P17252, P05771, Q05513, Q05655, Q04759, P05129, P05696, P24723
Entrez	5578, 18750
Expiration Date	12 months from date of receipt.