

Product Datasheet

Actin-pan Blocking Peptide (orb706241)

Description	The peptide is used to block Anti-Actin-pan antibody
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	ACTA2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62736, P62738, P68032, Q562R1, P63261, P63267, P62737, P68133, P60709
Entrez	11475, 81633, 59
Expiration Date	12 months from date of receipt.