

Product Datasheet

CCR5 Blocking Peptide (orb706297)

Description	The peptide is used to block Anti-CCR5 antibody <a href="https://www.biorbyt.com/ccr5-antibody-orb51..."
Reactivity	Human, Mouse, Primate, Rabbit
Tested Applications	Blocking
Target	CCR5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51682, P51681
Entrez	12774, 1234
Expiration Date	12 months from date of receipt.