

Product Datasheet

DAP Kinase 2 Blocking Peptide (orb706306)

Description The peptide is used to block Anti-DAP Kinase 2 antibody <a href="https://www.biorbyt.com/dapk2-antib..."

Reactivity Human, Mouse

Tested Applications Blocking

Target DAPK2

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **Q8VDF3, Q9UIK4**

Entrez **23604, 13143**

Expiration Date 12 months from date of receipt.