

Product Datasheet

Collagen-pan Blocking Peptide (orb706355)

Description	The peptide is used to block Anti-Collagen-pan antibody
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	COL1A1;COL2A1;COL3A1;COL5A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05539, P28481, P02454, P02461, Q9JI03, P13941, P02452, P08121, P11087, O88207, P20908, P02458
Entrez	25412, 1280, 85490, 1277, 12825, 12824, 12831, 1289, 29393, 1281, 84032, 12842
Expiration Date	12 months from date of receipt.