

Product Datasheet

Histone H2B (AcK85) Blocking Peptide (orb706359)

Description

The peptide is used to block Anti-Histone

[href=https://www.biorbyt.com/histo...](https://www.biorbyt.com/histo...)

Reactivity

Human, Mouse, Rat

Tested

Blocking

Applications
Target

HIST1H2BC; HIST1H2BE; HIST1H2BF;
HIST1H2BG; HIST1H2BI; HIST1H2BD;
HIST1H2BH; HIST1H2BK; HIST1H2BL;
HIST1H2BM; HIST1H2BN; HIST2H2BF;
H2BFS

Form/Appearance

Lyophilized powder

Storage

Shipped at 4°C. Store at -20°C for one year.

Note

For research use only

Application notes

Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity

> 85%

Source

Synthetic

Uniprot ID

P58876, Q64478, P62807, P33778, Q6ZWY9, Q99877, P23527, P06899, O60814, Q93079, Q99879, Q64475, Q99880

Entrez

8346, 85236, 8339, 319182, 8343, 8347, 440689, 8342, 319179, 3017, 319178, 8341, 8340, 8345, 8344

Expiration Date

12 months from date of receipt.