

Product Datasheet

JNK1/2/3 (phospho-Y185) Blocking Peptide (orb706403)

Description The peptide is used to block Anti-JNK1/2/3 (pY185) antibody

Reactivity Human, Mouse, Rat

Tested Applications Blocking

Target MAPK8; MAPK9; MAPK10

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID [Q61831](#), [P45983](#), [P49186](#), [P49187](#), [P49185](#), [Q9WTU6](#), [P45984](#), [P53779](#), [Q91Y86](#)

Entrez [5602](#), [50658](#), [26419](#), [26420](#), [5601](#), [25272](#), [5599](#), [116554](#)

Expiration Date 12 months from date of receipt.