

Product Datasheet

JNK1/2/3 Blocking Peptide (orb706422)

Description	The peptide is used to block Anti-JNK1/2/3 antibody
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAPK8; MAPK9; MAPK10
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61831, P45983, P49186, P49187, P49185, Q9WTU6, P45984, P53779, Q91Y86
Entrez	5602, 50658, 26419, 26420, 5601, 25272, 5599, 116554
Expiration Date	12 months from date of receipt.